

Bio Bums

High Absorbent Inserts

Australian made Bio Bums High Absorbent Inserts are comprised from the highest quality all natural & certified organic fabrics, with a unique folding system aimed to accommodate all types of wetters including boys, girls, heavy wetters & newborns. Two layers of highly absorbent silky soft Bamboo & Organic Cotton fabric with a thick and thirsty three layered fold out Hemp booster.

Unique Two Stage Folding System

Available in sizes Small & Large for use with newborn babies up to toddlers weighing 16kg and over. The unique Bio Bums Two Stage Folding System allows the insert to adjust in size accordingly as your baby grows which provides the longest nappy life from your High Absorbent Inserts.

Ideal Uses

Compatible with a variety of nappy systems including pocket nappies, fitted & pouch nappy covers, AIO & A12 Nappies & toddler toilet training pants.

Key Features

- ✓ High Absorbent
- ✓ Natural & Organic Fabrics
- ✓ Quick Dry Design
- ✓ Unique Folding System
 - ✓ Easy to use
 - ✓ Economical
- ✓ Environmentally Friendly
- ✓ Australian Made & Owned



Australian Made
Economical & Eco Friendly



Bio Bums offer families a fabulous range of nappy products that are 100% Australian made.

The fabrics selected and used in our products are made from natural and organic ingredients where ever possible, whilst having minimal impact on the environment.

Bio Bums inserts & Boosters are superior in absorbency and design. In addition to the high absorbent fabrics used, the patent pending design makes them quicker drying and easier to use.


When parents are deciding on which nappy system to use for their baby, Bio Bums inserts make the perfect choice for absorbency as they are compatible to use with most brands of nappy covers & nappy systems.

Bio Bums

P.O. Box 8338
Tarneit, Victoria 3029

Wholesale Enquiries: 0439 711 785
Phone: 03 8683 9579 or 03 8683 9079
E-mail: info@biobums.com.au
www.biobums.com.au

Bio Bums




The
Natural Nappy
Choice


High Absorbent Inserts
Product
Information

TRADITIONAL FOLD


- 1 Open & lay insert onto a flat surface.



- 2 Fold the hemp booster upwards onto the bamboo section of the insert.



- 3 Fold the right 1/3 of the bamboo towards the middle & over the hemp booster.




- 4 Fold the left side 1/3 of the bamboo body over the already pre-folded booster section.




High Absorbent Insert is now ready to be used with one of the many compatible nappy covers.

BOY FOLD


- 1 Open & lay insert onto a flat surface.




- 2 Rotate in an anticlockwise direction so that the hemp booster is positioned to the right of the bamboo body.



- 3 Fold the hemp booster onto the middle third of the bamboo part of the insert.



- 4 Fold the hemp booster at a 45° angle starting from the top corner where it is attached to the bamboo. Once folded, ensure the booster edge sits parallel with the bamboo section of the insert.




- 5 Fold the left 1/2 of the bamboo body over the folded booster section.


- 6 Position the insert with the folded hemp booster to the front of a nappy cover. Insert is now ready to use for baby boy, having concentrated absorbency towards the front where it is needed most.

GIRL FOLD



- 1 Open & lay insert onto a flat surface.




- 2 Fold the hemp booster upwards onto the bamboo body of the insert.



- 3 Fold the hemp booster in half so that it is doubled in thickness.



- 4 Rotate the insert in an anticlockwise direction so that the hemp booster is positioned to the right side of the bamboo body.



- 5 Fold the left 1/2 of the bamboo body over the folded booster section.

- 6 Place the insert in a nappy cover. This fold provides greater absorbency in the mid section of a nappy which is ideal to use for a baby girl.

Bio Bums

TWO STAGE FOLDING SYSTEM

Small High Absorbent Inserts - 450ml Liquid Capacity

Stage One Traditional Fold - Newborns & babies up to 6kgs

Stage Two Boy or Girl Fold - Babies from 6 - 16kgs

Large High Absorbent Inserts - 600ml Liquid Capacity

Stage One Traditional Fold - For babies 6kgs and above